

# Steps

1) Bind

## Plasmid DNA Prep

- 1) Resuspend bacteria in 250 uL P1
- 2) Lyse with 250 uL P2, wait 1-5 min
- 3) Neutralize with 350 uL chilled N3
- 4) Spin 5 min, load 600-650 uL SN on column
- 5) Incubate at RT for 2 min
- 5) Spin, discard FT

## Gel Extraction

- 1) Dissolve every 100mg of gel slice with 300 uL QG  
**(No heat needed w/ low melt agarose, invert or vortex)  
(Some protocols suggest 400 uL/100mg, also works)**
- 2) Load column with 600 uL, incubate at RT for 2 min
- 3) Spin, discard FT
- 4) Repeat until all dissolved gel passed through column  
**Only first incubation at RT is necessary to activate silica**

## PCR/Enzymatic Reaction Purification

- 1) Add 5 volumes of PB to 1 volume PCR/Enz reaction
- 2) Load 600 uL onto column, incubate at RT for 2 min
- 3) Spin, discard FT
- 4) Repeat until all product passed through column  
**Only first incubation at RT is necessary to activate silica**

2) Wash

- 1) Add 600-700 uL 1X PE (**diluted with EtOH!!!**), spin, discard FT
- 2) Repeat for second wash
- 3) Discard final FT, spin additional 2 minutes
- 4) Place column into fresh 1.5 mL tube

3) Elute

- 1) Add 15-100uL of elution buffer, pre-warm to 50C for DNA >5-10 kb
- 2) Spin, check concentration, transfer to fresh tube if feeling fancy.  
**\*\*\* Elution with volumes >50 uL will likely require concentration, Centrivap is king here**

# Recipes

### Buffer N3

4M Guanidine-HCl  
0.5M Potassium Acetate pH 4.2

### Buffer QG

5.5 M Guanidine Thiocyanate  
20 mM Tris HCl pH 6.6

### Buffer PB

5 M Guanidine-HCl  
20 mM Tris-HCl pH 6.6  
30% ethanol

\*\*\* Refer to Plasmid Prep Cheat sheet for all recipes \*\*\*

1) Bind

2) Wash

3) Elute

### 5X Buffer PE

80 mM NaCl  
8mM Tris-HCl pH 7.5  
**DILUTE WITH ETHANOL**

### 1X Buffer PE

1 vol. 5X Buffer PE  
**4 vol. 100% Ethanol**

**Elution Buffer:** 5 mM Tris-HCl pH 8

V1.2 Nov 24/19

**SN** = Supernatant **W/** = with  
**FT** = Flow through **Vol.** = Volume  
**Spin** = Centrifuge at 16kxg for 30s