

Gel Purification Protocol

- 1) Pool PCR products from multiple reactions
 - Anywhere from 4-12 reactions
 - “Can I gel purify 1 reaction?” Sure, but expect low yield...good luck downstream!
- 2) Add 0.5-1uL of DpnI to pooled reaction per 100 uL of PCR product -> 37°C for 1+ hours (Overnight is fine)
- 3) Desalt/concentrate digested PCR product with PCR Cleanup Protocol (Silica column cheat sheet)
 - Just concentrate PCR product and run on gel? -> Smear on agarose due to concentrated salt/buffer
 - Large volume difficult to load on low-melt agarose gel -> Diluted in gel volume
 - Silica column desalts AND concentrates pooled PCR product (Better than G50 resin, trust)
- 4) Elute in 35-50 uL if you don't have a Centrivap, elute in 100 uL and concentrate to 35 uL if you do
- 5) Load onto low-melt agarose with a blue-light reactive dye
 - Ethidium bromide + UV -> Damages DNA -> Longer DNA fragments = More damage/time
 - Don't have a gel comb that will fit 35-50uL? Tape a few lanes together
- 6) Cut out desired band under blue light and purify with silica column (Silica column cheat sheet)
- 7) Elute in 35-50 uL if you don't have a Centrivap, elute in 100 uL and concentrate to 35-50 uL if you do
- 8) Determine concentration of fragments with spectrophotometer:
 - **Less than 10ng/uL?** Dickey, but doable. Consider additional PCR reactions to bulk up concentration
 - **10-30ng/uL?** Good enough depending on size of fragments
 - **30-100ng/uL?** Comfortable amount for most applications
 - **>100ng/uL?** You're laughing! Enough for repeats!



MVP <3 I luv u