

## Gibson/Goldengate 2+ Fragment prep

1) Prepare PCR with:

- Highly accurate polymerase (Q5 or equivalent)
- Long fragments? GC Enhancer Buffer VERY useful
- Shorter fragments? <2-3kb? May or may not need GC enhancer -> Lowers accuracy (a bit)
- LOW template amount (0.1 to 0.5 ng/25 uL PCR reaction) -> Make it easier for the Dpnl

2) Thermocycling settings:

**Initial Denaturation:**

95°C for 90 seconds

**30-40 Cycles of:**

95°C for 30 seconds

50–67°C for 30 seconds →

72°C for 30 seconds/kb

**Have a gradient PCR machine?**

4-12 reactions from 55-67°C

**No gradient PCR?**

Start at 59°C and pray

**Final Extension:**

72°C for 2X Cycling Extension Time

END



**LET THE PCR MACHINE COAST TO ROOM TEMP!**

**PCR products can survive for DAYS at RT**

**Your thermocycler is not a refrigerator!**

# Gibson/Goldengate 2+ Fragment QC

- 1) Run 5-8uL of PCR product on a 1% agarose
- 2) Do I have a product of expected size? (\*)
- 3) Is the PCR product pretty clean?



# Gibson Assembly

- Determine concentration of all fragments with spectrophotometer
- Calculate fmol/uL for all fragments (NEB Biocalculator -> dsDNA Mass to Moles)
- 1 vector fragment, 1 insert? 1:3 molar ratio
- 1 vector fragment, 1 small insert (100-500bp)? 1:5 to 1:10 ratio
- 1 vector, 2 or more fragments? 1:1:1 etc ratio
- Shoot for 200fmol of all fragments combined, 50fmol is grim but doable, 400+ fmol for 3+ fragment assemblies
- Which mix to use?

## 1-4 Fragments?

### 1-3 Fragments?

#### **5X TEDA Cloning Mix (1 mL)**

500 mM Tris-HCl pH 7.5 (0.5 mL of 1M stock)  
50 mM MgCl<sub>2</sub> (50 uL of 1M stock)  
50 mM dithiothreitol (100 uL of 0.5M stock)  
0.25 g of PEG 8000 (May need gentle heat)  
1 µl of 10 U/µl T5 exonuclease (NEB)



**Combine fragments +  
mastermix in 10-20uL volume**



**30C for 40 minutes**

#### **2X Gibson Assembly Master Mix**

405µl Isothermal Start Mix  
25µl 1M DTT  
20µl 25mM dNTPs  
50µl NAD<sup>+</sup> (NEB Cat. B9007S)  
1µl T5 exonuclease (NEB Cat. M0363S)  
31.25 µl Phusion High Fidelity DNA Polymerase (NEB Cat. M0530S)  
250 µl Taq Ligase (NEB Cat. M0208S)  
467.75 µl H<sub>2</sub>O  
Mix by pipetting gently. Make 100 µl aliquots.

#### **Isothermal Start Mix**

1.5g PEG8000  
3ml 1M Tris-HCl, pH 8.0  
150µl 2M MgCl<sub>2</sub>  
Put on tube rotator until PEG is in solution

Recipe/Optimization by Ethan Ford -> [Ethanomics blog](#)



**Combine fragments + mastermix in 10-20uL volume**



**50C for 60 minutes**

## 1-4 Fragments?

**Commercial  
Gibson Mix  
(NEB or Other  
brands)**



**Combine fragments +  
mastermix in 10-20uL volume**



**50C for 60 minutes**

## 1-5 Fragments?

**Commercial  
Hi-fi Mix (NEB)**



**Combine fragments +  
mastermix in 10-20uL volume**



**50C for 60 minutes**

- Mix chosen depends on number of fragments to assemble and budget. Gibson/Hifi mastermix has a notoriously short shelf life! Even at -80C